

## OXYGEN ATOM BEAM SOURCE OBS



OBS 40 on DN40 (O.D. 2.75") CF-flange

The Oxygen Atom Beam Source OBS is a thermal gas cracker which produces an ion-free oxygen gas beam. The OBS does not cause ion induced damage on the substrate. It exhibits a very compact design and is easy to install and to operate.


The OBS has been developed in collaboration with Dr. Karl G. Tschersich from the Institute of Thin Films and Interfaces (ISG), Research Centre Jülich (FZJ). It is an extensively tested and well characterized Oxygen source.

The heart of the OBS is a long cracking tube heated by a surrounding filament. Operation temperatures up to 1900°C provide an efficient thermal cracking of O<sub>2</sub> molecules within the tube. The high purity cracking tube is the only part of the OBS with direct contact to the oxygen gas and forms a narrow angle distributed gas beam ejected towards the specimen.


The narrow beam of the OBS provides high atomic oxygen flux rates at the sample position while keeping the O<sub>2</sub> background pressure of the chamber limited as compared to plasma sources. In addition, the high cracking efficiency of up to 80%, which is available for low flux applications, reduces the chamber background pressure significantly as compared to ozone generators, where the ozone concentration is limited to about 10-15%.

Despite the high temperatures used for thermal gas cracking, the thermal load of the chamber is low due to the integrated water cooling and the use of a small diameter aperture.

- O<sub>2</sub> dissociation up to more than 80% depending on operation conditions
- Atomic O-flux density up to 1\*10<sup>15</sup>/cm<sup>2</sup>\*s
- No high-energy particles and ions
- Low power consumption
- Integrated water cooling, low thermal load on other experimental equipment



Grade of dissociation as a function of the gas flow for typical cracking tube temperatures of the OBS. Values are calculated for thermo-dynamical equilibrium condition.



Total flux of O-atoms leaving the cell orifice for different cell operation parameters (gas flow and tube temperature). A flux of about 1\*10<sup>17</sup> atoms/s corresponds to a typical deposition rate of 1 monolayer/s at a typical substrate distance of about 200mm.

## Applications

The use of the cracking tube results in a narrow angle distribution of the atomic Oxygen-beam (typical FWHM: 10-30°) which makes the OBS ideally suited for medium and low gas flux applications (up to 0.1 sccm) and smaller sample sizes or long substrate distances.

Due to the highly efficient oxygen cracking mechanism and the resulting low gas load of the vacuum system the OBS is well suited for use in standard UHV, MBE and ALD systems as a reliable, highly efficient and low cost atomic oxygen source.


Typical applications for the OBS are:

- Oxide layer deposition
- Surface cleaning procedures (e.g. removal of surface carbon contaminations)
- Surface oxidation (e.g. in ALD technology to form well defined oxide surface layers)
- Atomic spectroscopy of single atoms

## Technical Data

<b>Mounting Flange</b>	DN40CF (O.D.2.75")
<b>Dimensions in Vacuum</b>	190-400mm
<b>Gas injection</b>	Cracking tube
<b>Thermocouple</b>	W5%Re/W26%Re (type C)
<b>Bakeout temperature</b>	300°C
<b>Operating temperature</b>	up to 1900°C
<b>Cooling</b>	integrated water cooling
<b>Options</b>	Aperture (A), integrated Shutter (S)

Schematic drawing of the Oxygen Atom Beam Source OBS (Drawing shows OBS 40 with aperture and shutter)



## References:

- [1] K.G. Tschersich and V. von Bonin, J. Appl. Phys. 84, 4065 (1998)  
 [2] K.G. Tschersich, J. Appl. Phys. 87, 2565 (2000)

Dr. Eberl MBE-Komponenten GmbH  
 Gutenbergstrasse 8  
 71263 Weil der Stadt, Germany  
 Phone : +49 7033 6937-0  
 Fax : +49 7033 6937-20



info@mbe-components.com  
 www.mbe-components.com

USA / Canada  
 United Mineral & Chemical Corporation  
 Lyndhurst NJ 07071 (USA)  
 Phone : 1 800 777 0505  
 Fax : 201 507 1506  
 www.umccorp.com/

India  
 MACK INTERNATIONAL  
 Colaba, Mumbai 400 005  
 Phone : 22 2285 5261  
 Fax : 22 2285 2326  
 darius@mack.in

China  
 GERMAN TECH CO. Ltd.  
 Beijing 100083  
 Phone : +86 10 82867920  
 Fax : +86 10 82867919  
 www.germantech.com.cn/

Italy  
 OMICRON NanoTechnology GmbH  
 00136 Roma  
 Phone : 06 35 45 85 53  
 Fax : 06 35 40 38 67  
 omicron.italia@tiscali.it

Japan / Korea  
 VIEETECH JAPAN CO., Ltd  
 Mito, Ibaraki 311-1112  
 Phone : 29 269 9211  
 Fax : 29 269 9814  
 www.vieetech.co.jp/

Singapore  
 RESEARCH INSTRUMENTS PTE LTD  
 Singapore 139944  
 Phone : 6775-7284  
 Fax : 6775-9228  
 www.ri.com.sg/